




Dr. Overtoom Easy1 Bladder Access and Exchange procedure pack

Indication: *Difficultly accessible bladders*

- All the required materials in one single pack
- Reduced OR time
- No need for fluoroscopic verification due to use of a unique patented position verification catheter
- Unique guidewire guidable Foley catheter for ease of insertion



MEDpro Medical b.v.

Office Address:
Landjuweel 11
NL-3905 PE Veenendaal
The Netherlands

[T] +31 318 769 080
[F] +31 318 769 107
[E] info@medpromedical.nl

Chamber of commerce no.: 51758644

MKT/BRO/UR-0001-R0

www.medpromedical.nl



Dr. Overtoom Easy1 Bladder Access and Exchange procedure pack

Ordering information

Part number: PP-UR-0001

Kit contents:

- 1 x 14 Fr Foley catheter 35 cm
- 1 x Fascial Dilator 10 Fr
- 1 x Fascial Dilator 14 Fr
- 1 x Nitinol Guidewire hydro 0.035 / 80 cm
- 1 x Syringe for lubricating guide wire
- 1 x Syringe with 6 mL catheter lubricating gel
- 1 x Syringe with 10 mL balloon inflating gel
- 1 x Unique patented position verification catheter

MEDpro Medical b.v.

Office Address:
Landjuweel 11
NL-3905 PE Veenendaal
The Netherlands

[T] +31 318 769 080
[F] +31 318 769 107
[E] info@medpromedical.nl

Chamber of commerce no.: 51758644

MKT/BRO/UR-0001-R0

www.medpromedical.nl